

# Endo Pilot 2



Non-stop to a perfect Endo



EndoPilot<sup>2</sup> Set "premium"  
artitel-No.: 110 0106 .....



## The successful complete solution - EndoPilot - newly designed

As a modern company, it was our aim to present our successful, highly valued EndoPilot series in a fresh design. Look forward to tried and tested, sophisticated technology with a clear menu for your root canal treatments with the **EndoPilot<sup>2</sup>**.

The EndoPilot<sup>2</sup> range offers proven, high-quality, safe and efficient root canal treatment, from preparation to final thermoplastic filling according to the latest endodontic requirements. All necessary steps of a root canal treatment can be carried out with **only one** comfortable, space-saving control unit.

The EndoPilot<sup>2</sup> series has also been developed modularly, it is available and upgradeable in various configurations. Up to **7 functions** are possible. An upgrade, for e.g. with Obturation can be easily done by the user.



all 1300 38 66 22  
or ree Demo To a

- Modern, stylish design
- 7" color touch screen with front glass
- Large lithium-ion battery (up to 16 hours operational time)
- Fast charging function
- Easy software update via microSD card
- DownPack connection fully integrated (no changeover required)
- Improved sound output with sound effects
- Concealed cable routing
- High quality metal holder
- Bluetooth wireless footswitch
- New intuitive menu navigation
- Modularly expandable
- Flexible to use
- Compatible with predecessor model



Apex




Motor




DownPack



BackFill



UltraSonic



Battery



Wireless



Pump



Exclusive Australian Distributor  
T: 1300 38 66 22  
E: [enquiries@eurobiotechnologies.com.au](mailto:enquiries@eurobiotechnologies.com.au)  
W: [www.eurobiotechnologies.com.au](http://www.eurobiotechnologies.com.au)



Schlumbohm 50 Years  
1968-2018



### EndoPilot² Set „comfort“

Article no.: 110 0101 .....

Apex measurement + EndoMotor

**Control unit:**.....

**EndoPilot²**, Touch-Display,  
with Li-ion high performance battery  
Med. **AC/DC converter** 100-240V/50-60Hz  
**Apex cable** with File-clamp and Retainer  
**Wireless footswitch**, single

**Preparation :** .....

**EDP-Motor** ISO-E coupling, LED-indicator,  
integrated apex-connection  
**Contra-angle**, 1:1, for apex measurement,  
fully insulated, ISO-E coupling  
**Motor holder**, right

**Upgradeable:**

- DownPack
- BackFill
- UltraSonic
- Peristaltic-Pump
- SAF-System



### EndoPilot² Set „plus“

Article no.: 110 0102 .....

Apex measurement + EndoMotor +  
Obturation

**Control unit:**.....

**EndoPilot²**, Touch-Display  
with Li-ion high performance battery  
Med. **AC/DC converter** 100-240V/50-60Hz  
**Apex cable** with File-clamp and Retainer  
**Wireless footswitch**, single

**Preparation :** .....

**EDP-Motor** ISO-E coupling, LED-indicator,  
integrated apex-connection  
**Contra-angle**, 1:1, for apex measurement,  
fully insulated, ISO-E coupling  
**Motor holder**, right

**Obturation:** .....

**DownPack** handpiece, with LED-indicator  
**Heating tip** size: fine  
**Back-Fill Gun Obtura**  
**Gutta Percha Bar** (100 pcs/box)  
**Heat Insulator** (4 pcs/box)  
**Needles** Ø25 ga ( 5 pcs/box)  
**Accessory tools and cleaning set**  
**Handpiece holder**, left

**Upgradeable:**

- UltraSonic
- Peristaltic-Pump
- SAF-System



### EndoPilot² Set „ultra“

Article no.: 110 0103 .....

Apex measurement + EndoMotor +  
UltraSonic

**Control unit:**.....

**EndoPilot²**, Touch-Display  
with Li-ion high performance battery  
Med. **AC/DC converter** 100-240V/50-60Hz  
**Apex cable** with File-clamp and Retainer  
**Wireless footswitch**, single

**Preparation:** .....

**EDP-Motor** ISO-E coupling, LED-indicator,  
integrated apex-connection  
**Contra-angle**, 1:1, for apex measurement,  
fully insulated, ISO-E coupling

**UltraSonic:** .....

**Ultrasonic modul** with twin holder  
**Ultrasonic handpiece** Satelec Newtron  
with instrument tool

**Upgradeable:**

- Obturation
- Peristaltic-Pump
- SAF-System



### EndoPilot² Set „ultra plus“

Article no.: 110 0104 .....

Apex measurement + EndoMotor +  
UltraSonic + Obturation

**Control unit:**.....

**EndoPilot²**, Touch-Display  
with Li-ion high performance battery  
Med. **AC/DC converter** 100-240V/50-60Hz  
**Apex cable** with File-clamp and Retainer  
**Wireless footswitch**, single

**Preparation:** .....

**EDP-Motor** ISO-E coupling, LED-indicator,  
integrated apex-connection  
**Contra-angle**, 1:1, for apex measurement,  
fully insulated, ISO-E coupling

**UltraSonic:** .....

**Ultrasonic modul** with twin holder  
**Ultrasonic handpiece** Satelec Newtron  
with instrument tool

**Obturation:** .....

**DownPack** handpiece, with LED-indicator  
**Heating tip** size: fine  
**Back-Fill Gun Obtura**  
**Gutta Percha Bar** (100 pcs/box)  
**Heat Insulator** (4 pcs/box)  
**Needles** Ø25 ga ( 5 pcs/box)  
**Accessory tools and cleaning set**  
**Handpiece holder**, left

**Upgradeable:**

- Peristaltic-Pump
- SAF-System